


# Power Mic II Buttons (front)



**'Transfer Text' button:** clicking this allows you to transfer text from the dictation box into a free-text area of your EMR. Make sure your mouse cursor is clicked in a free-text area prior to clicking the 'Transfer text' button. You may also say 'Transfer Text'.

**'Previous field' button:** clicking this allows one to navigate to the previous variable bracket item when using a templates with brackets. Note: Saying 'previous field' does the same thing.

**Rewind button:** Not currently activated.

**Stop/Play button:** Puts the dictation box on/off Hidden mode. Works as a toggle. Allows one to hide Dictation Box while dictating, while reviewing areas of the chart. (e.g. labs, radiology reports, etc.) Can bring Dictation box back into view by clicking same button.

**Custom Left button:** Not currently activated.

**Left-click button:** Functions like the left-click on a mouse.

**Press to Talk button:** Keep this button pressed while dictating.

**'Next field' button:** clicking this allows one to navigate to the next bracket item, when using templates with brackets. Note: You can also say 'Next field' to do the same thing.

**Fast-forward button:** Opens the Dictation Box. Can also say "Show Dictation Box" to open it.

**'Accept Defaults' button:** This button removes all remaining brackets from the text (when using templates with brackets).

**Custom Right button:** Not currently activated.

**Right-click button:** Functions like the right-click on a mouse.

**Pointer:** Functions like a mouse pointer when rolled.

**Bottom Radio button:** Not currently activated.